



NICO® Fluid System

Simple solution to irrigating the surgical field

The NICO Fluid System is part of the NICO family of products that provide surgeons and OR staff with simple solutions to unmet needs in the surgical environment.

The NICO Fluid System provides controlled irrigation to the surgical site and refills automatically after each syringe compression. The extended tubing allows for fluid to be directly drawn from the warming bath to deliver fluid at a constant temperature.

- Delivers fluid to the surgical field
- Refills automatically with fluid after each syringe compression
- Eliminates the need to constantly exchange syringes



The NICO Fluid System is a syringe that comes pre-assembled and attached to a tubing set that includes a spike for use with saline or other fluid bag.

The NICO Fluid System is Product #NN-8015
To order, call 888.632.7071



The NICO Fluid System is not intended to replace the bulb syringe, typically used for gross flushing/irrigation of debris from the surgical site.



LIT1831RevC